# **CDN-ME-1902 – Multi-Element**

## SUMMARY

The application note summarizes the digestion of CDN-ME-1902, a multi-element Reference Material using ColdBlock<sup>™</sup> Digestion Pro Series Technology.

Instrument:	ColdBlock CBM sample digester, chiller, HF compatible liners, ICP-MS & ICP-OES				
Published:	January 2023				
Digestion Time:	30 Minutes				
Acid Used:	Reverse Aqua Regia, HF & $H_3BO_3$				
Average ColdBlock Recovery vs. CRM:	<ul> <li>102% Copper</li> <li>104% Lead</li> <li>103% Zinc</li> </ul>				

## METHODOLOGY

- 1. Chiller temperature was set to -5°C
- 2. 0.25g of each sample was weighed and placed into a ColdBlock™ Digestion vessel
- 3. 20 mL of reverse Aqua Regia + 3 mL HF was added
- 4. Sample was digested at 80% power for 20 minutes
- 5. 20mL of 4% , Boric acid was added
- 6. Samples were digested again at 80% power for 10 minutes
- Samples were cooled and bulked to 50mL using 2% HNO<sub>3</sub> + 0.5% HCl<sub>v/v</sub>

## DISCUSSION

- The addition of Boric acid will help re-solubilize any insoluble fluorides and will help neutralize any remaining HF in solution
- To improve silver recoveries, bulk up with >20% HCl<sub>v/v</sub>
- If the Sulfide content of your sample is > 10 wt.% reverse the ratios of Aqua Regia and use 1:3, HCI:HNO<sub>3</sub> - always add the Nitric acid first (reddish brown NO<sub>2</sub> fumes might form)



Standard CDN-ME-1902 was prepared by combining miscellaneous ores. cdnlabs.com

CDN-ME-1902; Multi-Element; CDN Resource Laboratories Ltd; Langley, British Columbia (November, 2019)

## CDN-ME-1902 - Multi-Element

#### Results

CDN-ME-1902												
Method:	0.25g	20mL reverse Aqua Regia + 3 mL HF digested at 80% for 20 minutes, add 20mL of 4% Boric Acid - and digest again at 80% for another 10 minutes										
Element	CDN Labs Certified 4-acid Value (ppm)	CDN Labs 95% Confidence Limits		Sample	Sample	Sample	Average	<u>cul</u>	%	% Recovery		
		Low	High	A	В	c	(ppm)	Stdev	RSD	vs 4-acid value		
Ag	349	332	366	301	301	310	304	4.06	1.3%	87%		
Cu	7810	7540	8080	7949	7986	8004	7980	23.09	0.3%	102%		
Pb	22000	21000	23000	22648	22768	23041	22819	164.59	0.7%	104%		
Zn	36600	34300	38900	37396	37827	37738	37654	185.77	0.5%	103%		

